

Product Information

VACUETTE® Glucose Tubes

General Information

Application

VACUETTE® Glucose Tubes are used for determinations in stabilised whole blood or plasma for glucose testing.

Materials

VACUETTE® Glucose Tubes are made out of PET (polyethyleneterephthalate). The caps are produced from PE (polyethylene) and are coloured by pigments purified in accordance with EN 71/3. The rubber component of the cap is composed of Brom Butyl Caoutchouc. The stabilisation ring is composed of PP (polypropylene).

Manufactured by

Greiner Bio-One GmbH, Austria

Certified according to SN EN ISO 9001:2008 and SN EN ISO 13485:2003 + AC 2007

Sterility

Sterile interior: SAL 10^{-6} (SAL = Sterility Assurance Level)

Standards: EN 552, EN 556, ANSI/AAMI/ISO 11137

Conformity

ISO 6710 "Single-use containers for venous blood specimen collection"

H1-A5 "Evacuated Tubes and Additives for Blood Specimen Collection- 5th Edition"; Approved Standard

CE mark

IVDD 98/79/EC "Directive 98/79/EC of the European Parliament and council of 27th of October 1998 on in vitro diagnostic medical devices" Classification: Other device (all devices except Annex II and self-testing devices).

Packaging

Tubes are packed in 50 pieces per rack wrapped with a polyethylene foil. Tubes are packed in 24 racks of 50 tubes. Each box is filled with 1200 tubes.


Shelf life

18 months

Product Information

VACUETTE® Glucose Tubes

General Information



| Labelling | Tube Label | Rack Label | Outer Carton |
|-----------------------|------------|------------|--------------|
| Manufacturers logo | X | X | X |
| VACUETTE® logo | X | X | X |
| Item number | X | X | X |
| Lot number | X | X | X |
| Expiry date | X | X | X |
| Fill volume | X | X | X |
| Fill line | X | | |
| Tube dimension | | X | |
| CE mark | X | X | X |
| Sterile R mark | X | X | X |
| Single use sign | X | X | X |
| Sterility indicator | | | X |
| Additive description | X | X | X |
| Cap description | | | X |
| Packaging information | | X | X |
| Storage guidelines | | X | X |
| Barcode label | (X) | X | X |

Wear gloves during venipuncture and when handling blood collection tubes to minimise exposure hazard.